

Technical Data

CV Reverse Osmosis System

The revolutionary R/O system



| | |
|---------------------------|--|
| Model | CV-RO-1000 |
| Order number | 8500499 |
| Time to market | 2016 |
| Description | Reverse Osmosis System |
| Certification | CE |
| Voltage | 100-240V |
| Frequency | 50,60Hz |
| Wattage | 45W |
| Environment temperature | 5 °C – 40 °C |
| Environment humidity | ≤ 80% |
| Water input quality | Tap water |
| Water input pressure | > 0.1Mpa |
| Water input temperature | 5 °C – 40 °C |
| Filters | Sediment filter - replace every 6 months Precarbon filter - replace every 6 months R/O membrane filter - replace every 30 months |
| Ability | 8 liters / hour |
| Net weight | 17kg |
| Gross weight | 20kg |
| Product dimensions | 440mm x 400mm x 140mm |
| Carton dimensions | 450mm x 450mm x 150mm |
| Installation | RO system can be installed wall mount or placed on table in horizontal or vertical position |
| Quality of purified water | ≤20µS/cm |
| Screen | Alphanumeric display |
| Packing list | 1x RO water purifier, 1x Sediment filter, 1x Precarbon filter, 1x R/O membrane filter, 1x Fixing plate for wall, 3x water pipe, 1x gooseneck tap, 1x power cable |
| Manual languages | English, German |
| Warranty | 2 years |

Certoclav Sterilizer GmbH
Georg Grinninger Str. 37
A-4050 Traun
Tel: +43 (0)7229 689
www.certoclav.com

Orders, technical assistance, custom orders, retailer contacts:
Christina Bokesch, Tel: +43 (0)7229 689 26, Fax: +43 (0)7229 689-20
Email: bokesch@certoclav.com
Product manager: Klaus David, Tel: +43 (0)7229 689 28, Email: david@certoclav.com
Manuals, technical information and instruction videos: www.certoclav.com