

Vertical laboratory autoclaves

CertoClav laboratory autoclave EL (Thermostat controlled)							
Ref.no.	Item	Programs	Temperature	Capacity	Chamber	Dimensions	Price €
8501312	CertoClav EL 12 GS 125/140 °C	2 fixed	125°C or 140°C	12 l	Aluminum	24 x 26 cm	2.040,00
8501412	CertoClav EL 12 GS 125/140 °C – SET with wire basket	2 fixed	125°C or 140°C	12 l	Aluminum	24 x 26 cm	2.148,00
8501312A	CertoClav EL 12 GS 115/121 °C	2 fixed	115°C or 121°C	12 l	Aluminum	24 x 26 cm	2.040,00
8501412A	CertoClav EL 12 GS 115/121 °C - SET with wire basket	2 fixed	115°C or 121°C	12 l	Aluminum	24 x 26 cm	2.148,00
8501318	CertoClav EL 18 GS 125/140 °C	2 fixed	125°C or 140°C	18 l	Aluminum	24 x 39 cm	2.467,00
8501418	CertoClav EL 18 GS 125/140 °C - SET with wire basket	2 fixed	125°C or 140°C	18 l	Aluminum	24 x 39 cm	2.540,00
8501318A	CertoClav EL 18 GS 115/121 °C	2 fixed	115°C or 121°C	18 l	Aluminum	24 x 39 cm	2.467,00
8501418A	CertoClav EL 18 GS 115/121 °C - SET with wire basket	2 fixed	115°C or 121°C	18 l	Aluminum	24 x 39 cm	2.540,00

CertoClav laboratory autoclave MultiControl (Microprocessor controlled with PT100 Sensor)							
Ref.no.	Item	Programs	Temperature	Capacity	Chamber	Dimensions	Price €
8502312A	CertoClav MultiControl 12 l	6 fixed 2 programmable	selectable (max. 140°C)	12 l	Aluminum	24 x 26 cm	3.577,00
8502412A	CertoClav MultiControl 12 l - SET with wire basket	6 fixed 2 programmable	selectable (max. 140°C)	12 l	Aluminum	24 x 26 cm	3.643,00
8502318A	CertoClav MultiControl 18 l	6 fixed 2 programmable	selectable (max. 140°C)	18 l	Aluminum	24 x 39 cm	3.928,00
8502418A	CertoClav MultiControl 18 l - SET with wire basket	6 fixed 2 programmable	selectable (max. 140°C)	18 l	Aluminum	24 x 39 cm	3.996,00

CertoClav laboratory autoclave Essential							
Ref.no.	Item						
8510173	CertoClav Essential 12 l	REPLACED with CertoClav MultiControl (Ref.no. 8502312A)				NOT AVAILABLE ANYMORE	
8510183	CertoClav Essential 18 l	REPLACED with CertoClav MultiControl (Ref.no. 8502318A)				NOT AVAILABLE ANYMORE	

Accessories for vertical laboratory autoclaves							
Ref.no.	Item					Accessories for	Price €
8583010	Wire basket with handle	ø 23 cm		Stainless steel 18/10		EL, MultiControl	158,90
1250499	Service Set 1	lid gasket, valve gasket, gasket for exhaust cock and manometer, steam valve		gaskets for a lid maintenance		EL, MultiControl	44,00
8582020	Condensation vessel	capacity: 5 l		with exhaust hose		EL, MultiControl	276,90
8582016	Exhaust hose	length: 1 m				Condensation vessel (8582020)	32,90
8500210	PC-Software Pro	allows multi-level programs (several time and temperature steps)		unlimited license key		MultiControl	242,90
1300187	Log printer	for logging during sterilization on thermal paper		external device		MultiControl	316,90
8500496	Thermal paper	10 rolls á 40 m				Log printer (1300187)	12,50
8582010	Exhaust filter	with condensation vessel and exhaust hose		0,01 Micron		MultiControl	500,50
8582011	Filter cartridge	0,01 Micron				replacement for exhaust filter (8582010)	264,50

Prices EXW excl. VAT

Deliveries are based to the terms and conditions published on our website.

Vertical laboratory autoclaves

CertoClav laboratory autoclave Classic (Thermostat controlled with timer)							
Ref.no.	Item	Programs	Temperature	Capacity	Chamber	Dimensions	Price €
8510174	CertoClav Classic 125/140 °C	2 fixed	125°C or 140°C	18 l	Stainless steel	26 x 35 cm	2.040,00
8510184	CertoClav Classic 115/121 °C	2 fixed	115°C or 121°C	18 l	Stainless steel	26 x 35 cm	2.040,00
CertoClav laboratory autoclave MultiControl 2 (freely programmable for any program, Microprocessor controlled with PT100 Sensor)							
Ref.no.	Item	Programs	Temperature	Capacity	Chamber	Dimensions	Price €
8510175	CertoClav MultiControl 2 with Controller MC2	6 fixed 2 programmable	selectable (max.140°C)	18 l	Stainless steel	26 x 35 cm	2.928,00
8510298	CertoClav MultiControl 2 Touch with Controller MC2, 7" Touchscreen, Cloud Basic	6 fixed 2 programmable	selectable (max.140°C)	18 l	Stainless steel	26 x 35 cm	3.288,00
CertoClav laboratory autoclave Connect (full automatic laboratory autoclav with cloud connection)							
Ref.no.	Item	Programs	Temperature	Capacity	Chamber	Dimensions	Price €
8510176	CertoClav Connect with 7" Touchscreen, Cloud Basic	6 fixed 2 programmable	selectable (max.140°C)	18 l	Stainless steel	25 x 35 cm	3.488,00
Accessories for vertical laboratory autoclaves							
Ref.no.	Item					Accessories for	Price €
1300149	Wire basket with handle	ø 25 cm		Stainless Steel 18/10		Classic, MultiControl 2, MultiControl 2 Touch, Connect	158,90
1310186	Service Set 2	lid gasket, control valve gasket, seals for exhaust cock, pressure gauge, steam valve, safety valve		gaskets for a lid maintenance		Classic, MultiControl 2, MultiControl 2 Touch, Connect	61,50
8582020	Condensation vessel	capacity: 5 l		with exhaust hose		Classic, MultiControl 2, MultiControl 2 Touch, Connect	276,90
8582016	Exhaust hose	length: 1 m				Condensation vessel (8582020)	32,90
8500210	PC-Software Pro	allows multi-level programs (several time and temperature steps)		unlimited license key		MultiControl 2	242,90
1300187	Log printer	for logging during sterilization on thermal paper		external device		MultiControl 2, MultiControl 2 Touch	316,90
8500496	Thermal paper	10 rolls á 40 m				Log printer (1300187)	12,50
8582010	Exhaust filter	with condensation vessel and exhaust hose		0,01 Micron		MultiControl 2, MultiControl 2 Touch	500,50
8582011	Filter cartridge	0,01 Micron				replacement for exhaust filter (8582010)	264,50
1300185	Media sensor	for the sterilization of liquids		only available in conjunction with a new device		Connect	363,90
8500508	CertoServer	local server with preinstalled autoclave management software for autoclave management and automatic backups				Connect	799,50
8500212	CertoCloud premium license	enables SMS notifications, storage for protocols, premium remote support, real-time-data online				Connect	263,90

Prices EXW excl. VAT

Deliveries are based to the terms and conditions published on our website.

Accessories for vertical autoclaves

Accessories for vertical autoclaves				
Ref.no.	Item		Accessories for	Price €
8583010	Wire basket with handle	ø 23 cm	Stainless Steel 18/10	Autoclaves aluminum 158,90
1300149	Wire basket with handle	ø 25 cm	Stainless Steel 18/10	Autoclaves S/S 158,90
1252203	Tray	ø 225 mm x 40 mm	Aluminum	Autoclaves aluminum 31,90
1300155	Tray	ø 225 mm x 40 mm	Stainless Steel	Autoclaves S/S 42,50
1250499	Service Set 1	lid gasket, valve gasket, gasket for exhaust cock and manometer, steam valve	gaskets for a lid maintenance	Autoclaves aluminum 44,00
1310186	Service Set 2	lid gasket, control valve gasket, seals for exhaust cock, pressure gauge, steam valve, safety valve	gaskets for a lid maintenance	Autoclaves S/S 61,50
8582020	Condensation vessel	capacity: 5 l	with exhaust hose	all vertical models 276,90
8582016	Exhaust hose	length: 1 m		Condensation vessel (8582020) 32,90
8500210	PC-Software Pro	allows multi-level programs (several time and temperature steps) unlimited license key		MultiControl 242,90
1300187	Log printer	for logging during sterilization on thermal paper	external device	MultiControl, MultiControl 2 316,90
8500496	Thermal paper	10 rolls á 40 m		Log printer (1300187) 12,50
8582010	Exhaust filter	with condensation vessel and exhaust hose	0,01 Micron	MultiControl, MultiControl 2 500,50
8582011	Filter cartridge	0,01 Micron		replacement for exhaust filter (8582010) 264,50
1300185	Media sensor	for the sterilization of liquids	only available in conjunction with a new device	Connect 363,90
8500508	CertoServer	local server with preinstalled autoclave management software for autoclave management and automatic backups		Connect 799,50
8500212	CertoCloud premium license	enables SMS notifications, storage for protocols, premium remote support, real-time-data online		Connect 263,90
1250000	Check and Service Package			all vertical models Net 160,00
1250001	Factory-made validation of the sterilization cycle		with protocol	all vertical models Net 260,00
8581000	Modification CertoClav EL	with new controller MC2/PT100		EL Net 1.640,00

Prices EXW excl. VAT

Deliveries are based to the terms and conditions published on our website.

Horizontal autoclaves

CertoClav laboratory autoclave Vacuum Pro								
Ref.no.	Item		Programs	Temperature	Capacity	Chamber	Dimensions	Price €
8500487	CertoClav Vac Pro 12	Stand-alone with internal water tanks	12 programmable	selectable (max. 134°C)	12 l	Stainless Steel	20 x 36 cm	3.793,00
8500487W	CertoClav Vac Pro 12	with fixed water connection	12 programmable	selectable (max. 134°C)	12 l	Stainless Steel	20 x 36 cm	3.793,00
8500488	CertoClav Vac Pro 18	Stand-alone with internal water tanks	12 programmable	selectable (max. 134°C)	18 l	Stainless Steel	25 x 35 cm	4.049,00
8500488W	CertoClav Vac Pro 18	with fixed water connection	12 programmable	selectable (max. 134°C)	18 l	Stainless Steel	25 x 35 cm	4.049,00
8500489	CertoClav Vac Pro 22	Stand-alone with internal water tanks	12 programmable	selectable (max. 134°C)	22 l	Stainless Steel	25 x 45 cm	4.305,00
8500489W	CertoClav Vac Pro 22	with fixed water connection	12 programmable	selectable (max. 134°C)	22 l	Stainless Steel	25 x 45 cm	4.305,00
8500610	CertoClav Vac Pro 45	Stand-alone with internal water tanks	8 programmable	selectable (max. 134°C)	45 l	Stainless Steel	31,9 x 61,7 cm	9.600,00

Prices EXW excl. VAT

Deliveries are based to the terms and conditions published on our website.

Accessories for horizontal laboratory autoclaves

Accessories for CertoClav Vacuum Pro				
Ref.no.	Item		Accessories for	Price €
8500508	CertoServer	local server with preinstalled autoclave management software for autoclave management and automatic backups	all models	799,50
8500212	CertoCloud premium license	enables SMS notifications, storage for protocols, premium remote support, real-time-data online	all models	263,90
8500513	Template for IQ/OQ documentation	on usb flash drive	Vac Pro 12-22	820,00
8500513	IQ/OQ qualification	Factory-made validation including documentation	Vac Pro 12-22	3.500,00
8500514	Exhaust filtration for 200 sterilization cycles	for safety sterilization of biological hazardous substances	only available in conjunction with a new device	all models 707,50

Prices EXW excl. VAT

Deliveries are based to the terms and conditions published on our website.

Accessories for horizontal laboratory autoclaves

Accessories for CertoClav Vacuum Pro					
Ref.no.	Item		Material:	Accessories for	Price €
8500479	Tray Set with 3 trays	middle tray 18 x 34 cm (W x L) upper and lower tray 14 x 34 cm (W x L)	Stainless steel	Vac Pro 12	205,00
8500480	Tray Set with 3 trays	middle tray 22 x 33 cm (W x L) upper and lower tray 17 x 33 cm (W x L)	Stainless steel	Vac Pro 18	266,50
8500481	Tray Set with 3 trays	middle tray 22 x 43 cm (W x L) upper and lower tray 17 x 43 cm (W x L)	Stainless steel	Vac Pro 22	276,90
8500404	Tray Set with 3 trays	middle tray 22 x 43 cm (W x L) upper and lower tray 17 x 43 cm (W x L)	Stainless steel	Vac Pro 45	315,90
8500479A	Tray for bottles	12,5 x 32 cm (W x L)	Stainless steel	Vac Pro 12	85,50
8500480A	Tray for bottles	16 x 30 cm (W x L)	Stainless steel	Vac Pro 18	95,90
8500481A	Tray for bottles	16 x 40 cm (W x L)	Stainless steel	Vac Pro 22	106,50
8500404A	Tray for bottles	16 x 40 cm (W x L)	Stainless steel	Vac Pro 45	114,50
8500401	Wire basket including drip tray	13,5 x 13,5 x 32 cm (W x H x L)	Stainless steel	Vac Pro 12	119,00
8500402	Wire basket including drip tray	17 x 17 x 31 cm (W x H x L)	Stainless steel	Vac Pro 18	159,00
8500403	Wire basket including drip tray	17 x 17 x 41 cm (W x H x L)	Stainless steel	Vac Pro 22	179,00
8500406	Wire basket including drip tray	17 x 17 x 41 cm (W x H x L)	Stainless steel	Vac Pro 45	376,00
8500495	Service Set 3	door sealing, bacteria filter, safety valve		Vac Pro 12-22	172,50
8500407	Service Set 4	door sealing, bacteria filter, safety valve		Vac Pro 45	197,90
1300187	Log printer	for logging during sterilization on thermal paper	external device	Vac Pro 12-22	316,90
8500496	Thermal paper	10 rolls á 40 m		Log printer (1300187)	12,50
8500482	Water distiller	distills fully automatic 1.5 liters of water within an hour	external device		215,00
8500499	CertoClav R/O-System including filter	reverse osmosis system for automatic water supply	requires tap water connection	models with fixed water connection	604,90
8500510	Sediment filter	replace every 6 months		CertoClav R/O-System	25,90
8500511	Precarbon filter	replace every 6 months		CertoClav R/O-System	25,90
8500512	Membrane filter	replace every 30 months		CertoClav R/O-System	102,50
8500502	Label printer	prints expiration date and protocol information on self-adhesive label	external device	Vac Pro 12-22	305,50
8500503	Barcode scanner	allows reading the sterilization protocol of an item	external device	Vac Pro 12-22	81,00
8500507	Direct thermal labels	10 rolls, self-adhesive, removable		Label printer (8500502)	19,50

Prices EXW excl. VAT

Deliveries are based to the terms and conditions published on our website.

Spare parts for **CertoClav** tabletop autoclaves

Gasket			
Ref.no.	Item	Spare Part for	Price €
1250404	loose Silicon white 24 cm for lid	Autoclaves aluminum	23,50
1260404	packed Silicon white 24 cm for lid	Autoclaves aluminum	24,00
1300122	white 26 cm for lid	Autoclaves S/S	43,90
1250433	for sterilization valve aluminum lid	Autoclaves aluminum	1,90
1300483	for sterilization valve stainless steel lid	Autoclaves S/S	2,50
1250437	for thermometer for CV-EL	CertoClav EL	3,90
1250451	for manometer and steam release cock	all vertical models	2,50
1250408	for base black 24 cm	Autoclaves aluminum	18,50
1300026	for base white 26 cm	Autoclaves S/S	24,50
Valves			
Ref.no.	Item	Spare Part for	Price €
1250304	CV-EL / Classic 125/140 °C	CertoClav EL CertoClav Classic 125/140 °C	121,00
1250302	CV-EL / Classic 115/121 °C	CertoClav EL CertoClav Classic 115/121 °C	121,00
1250320	for MultiControl / Essential 3,0 bar	CertoClav MultiControl CertoClav MultiControl 2 CertoClav Essential	121,00
1250242	Steam release valve aluminum lid	Autoclaves aluminum	30,50
1300123	Steam release valve stainless steel lid	Autoclaves S/S	32,90
1250310	Safety valve complete lid Aluminum	Autoclaves aluminum	30,90
1300103	Safety valve complete lid Stainless steel	Autoclaves S/S	33,90
Handles			
Ref.no.	Item	Spare Part for	Price €
1250235	Lid handle with lock complete	Autoclaves aluminum	44,50
1240239	Locking mechanism 4 parts only for repair	Autoclaves aluminum	31,50
1250236	Lid handle with screw and nut	Autoclaves aluminum	26,90
1250237	Handle with arrow for chamber Aluminum	Autoclaves aluminum	26,50
1250238	Handle without arrow for chamber Alu	Autoclaves aluminum	26,50
1250234	Lid handle with lock complete for ED	Autoclaves S/S	63,90
1300086	Lid handle with screw and nut for ED	Autoclaves S/S	22,90
1300012	Handle with arrow for chamber stainless-steel	Autoclaves S/S	19,50

Prices EXW excl. VAT

Deliveries are based to the terms and conditions published on our website.

Spare parts for CertoClav tabletop autoclaves

Lid			
Ref.no.	Item	Spare Part for	Price €
7940326	Lid Alu without parts CV-EL/MC	CertoClav EL CertoClav MultiControl	416,90
8511310	Lid CV-EL complete, without pressure control valve	CertoClav EL	1.034,50
8510312	Lid MultiControl, ES Alu complete with valve	CertoClav MultiControl CertoClav Essential	1.087,50
Parts			
Ref.no.	Item	Spare Part for	Price €
1250332	Thermometer with gasket	all vertical models	146,90
1250330	Glass scale for thermometer	all vertical models	59,90
1250334	Manometer EN837-1 Cl.1.6 with calibration	all vertical models	340,50
1250334A	Manometer EN837-1 Cl.1.6 with Manufact. Certificat acc. EN10204	all vertical models	211,50
1250348	Steam release cock CV-EL complete	CertoClav EL	126,90
1250194	Counter nut for manometer and steam release cock	all vertical models	2,90
1250208	Steam release cock MC with venting nozzle	CertoClav MultiControl CertoClav MultiControl 2	205,50
1240211	Plate with venting nozzle for MC 11 mm ø	CertoClav MultiControl CertoClav MultiControl 2	79,00
1250380	Cleaning needle for venting nozzle	all vertical models	24,90
Aluminum chamber with handles and base			
Ref.no.	Item	Spare Part for	Price €
8520112	CV-EL 12 l GS	CertoClav EL 12 l GS	1.420,50
8520118	CV-EL 18 l GS	CertoClav EL 18 l GS	1.927,50
8523312	MultiControl / Essential 12 l	CertoClav MultiControl 12 l CertoClav Essential 12 l	1.467,90
8523318	MultiControl / Essential 18 l	CertoClav MultiControl 18 l CertoClav Essential 18 l	1.974,50
Base with socket			
Ref.no.	Item	Spare Part for	Price €
1250324	CV-EL 12/18 l with or without fuses	CertoClav EL	126,50
1251324	MultiControl / Essential	CertoClav MultiControl CertoClav Essential	126,50
1300025	Classic / MC2	CertoClav Classic CertoClav MultiControl 2	143,50

Prices EXW excl. VAT

Deliveries are based to the terms and conditions published on our website.

Spare parts for CertoClav tabletop autoclaves

Cable set complete			
Ref.no.	Item	Spare Part for	Price €
1250357	for CV-EL without fuses	CertoClav EL	57,90
1250356	for MultiControl chamber Aluminum	CertoClav MultiControl	45,90
1310391	for Classic	CertoClav Classic	63,50
1310027	for MC2	CertoClav MultiControl 2	63,50
1250337	Cable with mains plug for CV-EL /Classic, 1250337UK= with UK-plug	all vertical models	48,50
Controller			
Ref.no.	Item	Spare Part for	Price €
1310243	Control unit MultiControl 2	CertoClav EL CertoClav MultiControl CertoClav MultiControl 2	1.254,30
1300188K	PT 100 Sensor with cable	CertoClav EL CertoClav MultiControl CertoClav MultiControl 2	226,50
1300197	Mains cable MultiControl 2/Essential	CertoClav EL CertoClav MultiControl CertoClav MultiControl 2	48,40
1300198	Heating cable MultiControl 2/Essential	CertoClav EL CertoClav MultiControl CertoClav MultiControl 2	48,40
1250360	Sensor cable green MC old version	CertoClav MultiControl	186,70
1250381	Mains cable 3 m MultiControl old version	CertoClav MultiControl	48,40
1250382	Heating cable 1,4 m MultiControl old version	CertoClav MultiControl	48,40

Prices EXW excl. VAT

Deliveries are based to the terms and conditions published on our website.

Spare parts for CertoClav tabletop autoclaves

Ref.no.	Item	Spare Part for	Price €
1250347	Green rocker switch with cover	Autoclaves aluminum	27,90
1250329	Tuning knob setting temp. CV-EL/Classic	CertoClav EL CertoClav Classic	13,50
1250338	Thermostat CV-EL Alu	CertoClav EL	56,90
1230338A	Thermostat Classic Stainless steel	CertoClav Classic	56,90
1250335	Bimetal fuse chamber Aluminum	Autoclaves aluminum	24,50
1300309	Bimetal fuse chamber stainless-steel	Autoclaves S/S	24,50
1250325	Socket complete for base black	Autoclaves aluminum	39,90
1240325A	Socket complete for base white	Autoclaves S/S	39,90
1250363	Fuse 10 Amp.	all vertical models	1,90
1300031	Timer Classic	CertoClav Classic	54,50
1300042	Knob for timer Classic	CertoClav Classic	50,50
1300036	Lamp for Classic	CertoClav Classic	51,90
Repairs			
Ref.no.	Item	Spare Part for	Price €
1300000	15 minutes working time	Use original CertoClav spare parts only for repairs. all autoclaves	Net 32,00

Prices EXW excl. VAT

Deliveries are based to the terms and conditions published on our website.



CERTOCLAV

Since 60 years.



Orders
Distributor verification
Custom-made products
Operating instructions and operating videos
Technical support



Certoclav Sterilizer GmbH
Peintner Straße 10
A-4060 Leonding



Phone: +43 732 674 278
E-Mail: support@certoclav.com
Website: www.certoclav.com