

**Vertical laboratory autoclaves**
**CertoClav laboratory autoclave EL**

Ref.no.	Item		Programs	Temperature	Capacity	Chamber	Dimensions	Price €
<b>8501312</b>	<b>CertoClav EL 12L</b>	Thermostat controlled	2 fixed	125°C or 140°C	12 l	Aluminum	24 x 26 cm	<b>1.999,00</b>
<b>8501312A</b>	<b>CertoClav EL 12L</b>	Thermostat controlled	2 fixed	115°C or 121°C	12 l	Aluminum	24 x 26 cm	<b>1.999,00</b>
<b>8501318</b>	<b>CertoClav EL 18L</b>	Thermostat controlled	2 fixed	125°C or 140°C	18 l	Aluminum	24 x 39 cm	<b>2.406,00</b>
<b>8501318A</b>	<b>CertoClav EL 18L</b>	Thermostat controlled	2 fixed	115°C or 121°C	18 l	Aluminum	24 x 39 cm	<b>2.406,00</b>
<b>8501412</b>	<b>CV EL 12L Set with wire basket</b>	Thermostat controlled	2 fixed	125°C or 140°C	12 l	Aluminum	24 x 26 cm	<b>2.095,00</b>
<b>8501412A</b>	<b>CV EL 12L Set with wire basket</b>	Thermostat controlled	2 fixed	115°C or 121°C	12 l	Aluminum	24 x 26 cm	<b>2.095,00</b>
<b>8501418</b>	<b>CV EL 18L Set with wire basket</b>	Thermostat controlled	2 fixed	125°C or 140°C	18 l	Aluminum	24 x 39 cm	<b>2.478,00</b>
<b>8501418A</b>	<b>CV EL 18L Set with wire basket</b>	Thermostat controlled	2 fixed	115°C or 121°C	18 l	Aluminum	24 x 39 cm	<b>2.478,00</b>

**CertoClav laboratory autoclave MultiControl**

Ref.no.	Item		Programs	Temperature	Capacity	Chamber	Dimensions	Price €
<b>8502312A</b>	<b>CertoClav MultiControl 12L</b>	Microprocessor controlled with PT100 Sensor	6 fixed 2 programmable	selectable (max. 140°C)	12 l	Aluminum	24 x 26 cm	<b>3.489,00</b>
<b>8502318A</b>	<b>CertoClav MultiControl 18L</b>	Microprocessor controlled with PT100 Sensor	6 fixed 2 programmable	selectable (max. 140°C)	18 l	Aluminum	24 x 39 cm	<b>3.832,00</b>
<b>8502412A</b>	<b>CV MultiControl 12L Set with wire basket</b>	Microprocessor controlled with PT100 Sensor	6 fixed 2 programmable	selectable (max. 140°C)	12 l	Aluminum	24 x 26 cm	<b>3.554,00</b>
<b>8502418A</b>	<b>CV MultiControl 18L Set with wire basket</b>	Microprocessor controlled with PT100 Sensor	6 fixed 2 programmable	selectable (max. 140°C)	18 l	Aluminum	24 x 39 cm	<b>3.898,00</b>

**CertoClav laboratory autoclave Essential**

Ref.no.	Item		Programs	Temperature	Capacity	Chamber	Dimensions	Price €
<b>8510173</b>	<b>CertoClav Essential 12L</b>	Microprocessor controlled with PT100 Sensor	1 fixed	121°C	12 l	Aluminum	24 x 26 cm	<b>2.078,00</b>
<b>8510183</b>	<b>CertoClav Essential 18L</b>	Microprocessor controlled with PT100 Sensor	1 fixed	121°C	18 l	Aluminum	24 x 39 cm	<b>2.250,00</b>

**Accessories for vertical laboratory autoclaves, Aluminum**

Ref.no.	Item				Accessories for	Price €
<b>8583010</b>	<b>wire basket with handle</b>	ø 23 cm		Stainless steel	EL, MC, Essential	<b>155,00</b>
<b>1250499</b>	<b>Service Set 1</b>	gaskets for a lid maintenance	with lid gasket, valve gasket, gasket for exhaust cock and manometer, steam valve		EL, MC, Essential	<b>42,60</b>
<b>1300187</b>	<b>log printer</b>	external	for logging during sterilization on thermal paper		MC, Essential	<b>309,00</b>
<b>8500496</b>	<b>thermal paper</b>	10 rolls á 40 m		suitable for log printer (1300187)	MC, Essential	<b>12,00</b>
<b>8582010</b>	<b>exhaust filter</b>	0,01 Micron		with condensation vessel and exhaust hose	MC, Essential	<b>488,00</b>
<b>8582011</b>	<b>filter cartridge</b>	0,01 Micron		replacement for exhaust filter (8582010)	MC, Essential	<b>258,00</b>
<b>8582016</b>	<b>exhaust hose</b>	length: 1 m		suitable for condensation vessel (8582020)	EL, MC, Essential	<b>32,00</b>
<b>8582020</b>	<b>condensation vessel</b>	capacity: 5 l		with exhaust hose	EL, MC, Essential	<b>270,00</b>
<b>8500210</b>	<b>PC software Pro license</b>		allows multi-level programs (several time and temperature steps)		MC	<b>237,00</b>

Prices EXW excl. VAT

Deliveries are based to the terms and conditions published on our website.

**Vertical laboratory autoclaves**
**CertoClav laboratory autoclave Classic**

Ref.no.	Item		Programs	Temperature	Capacity	Chamber	Dimensions	Price €
<b>8510174</b>	CertoClav Classic 18L	Thermostat controlled with timer	2 fixed	125°C or 140°C	18 l	Stainless steel	26 x 35 cm	<b>1.990,00</b>
<b>8510184</b>	CertoClav Classic 18L	Thermostat controlled with timer	2 fixed	115°C or 121°C	18 l	Stainless steel	26 x 35 cm	<b>1.990,00</b>

**CertoClav laboratory autoclave MultiControl 2**

Ref.no.	Item		Programs	Temperature	Capacity	Chamber	Dimensions	Price €
<b>8510175</b>	CertoClav MultiControl 2 18L	Microprocessor controlled with PT100 Sensor	6 fixed 2 programmable	selectable (max. 140°C)	18 l	Stainless steel	26 x 35 cm	<b>2.856,00</b>
<b>8510298</b>	CertoClav MultiControl 2 Touch 18L	Microprocessor controlled with PT100 Sensor	6 fixed 2 programmable	selectable (max. 140°C)	18 l	Stainless steel	26 x 35 cm	<b>3.207,00</b>

**CertoClav laboratory autoclave Connect**

Ref.no.	Item		Programs	Temperature	Capacity	Chamber	Dimensions	Price €
<b>8510176</b>	CertoClav Connect 18L	fully automatic 7" Touchscreen with cloud documentation	6 fixed 2 programmable	selectable (max. 140°C)	18 l	Stainless steel	25 x 35 cm	<b>3.402,00</b>

**CertoClav laboratory autoclave CultureMedia**

Ref.no.	Item		Programs	Temperature	Capacity	Chamber	Dimensions	Price €
<b>8510176</b>	CertoClav CultureMedia 18L	fully automatic integrated magnetic stirrer 7" Touchscreen	freely programmable	selectable (max. 140°C)	18 l	Stainless steel	26 x 35 cm	<b>7.500,00</b>

**Accessories for vertical laboratory autoclaves, Stainless steel**

Ref.no.	Item				Accessories for	Price €
<b>1300149</b>	wire basket with handle	ø 25 cm		Stainless steel	Classic, MC 2, MC 2 Touch, Connect, CM	<b>155,00</b>
<b>1310186</b>	Service Set 2	gaskets for a lid maintenance	with lid gasket, control valve gasket, seals for exhaust cock and pressure gauge, steam valve and safety valve		Classic, MC 2, MC 2 Touch, Connect, CM	<b>59,80</b>
<b>1300187</b>	log printer	External	for logging during sterilization on thermal paper		MC 2, MC 2 Touch	<b>309,00</b>
<b>8500496</b>	thermal paper	10 rolls á 40 m	suitable for log printer (1300187)		MC 2, MC 2 Touch	<b>12,00</b>
<b>8582010</b>	exhaust filter	0,01 Micron	with condensation vessel and exhaust hose		MC 2, MC 2 Touch, Connect, CM	<b>488,00</b>
<b>8582011</b>	filter cartridge	0,01 Micron	replacement for exhaust filter (8582010)		MC 2, MC 2 Touch, Connect, CM	<b>258,00</b>
<b>8582016</b>	exhaust hose	length: 1 m	suitable for condensation vessel (8582020)		Classic, MC 2, MC 2 Touch, Connect, CM	<b>32,00</b>
<b>8582020</b>	condensation vessel	capacity: 5 l	with exhaust hose		Classic, MC 2, MC 2 Touch, Connect, CM	<b>270,00</b>
<b>8500210</b>	PC software Pro license	allows multi-level programs (several time and temperature steps)			MC 2	<b>237,00</b>
<b>1300185</b>	media sensor		for the sterilization of liquids		Connect, CM	<b>355,00</b>
<b>8500508</b>	CertoServer		communicates with the autoclave via network automatic data backups & browser-based data access		Connect	<b>499,00</b>

Prices EXW excl. VAT

Deliveries are based to the terms and conditions published on our website.

## Horizontal laboratory autoclaves

CertoClav laboratory autoclave Vacuum Pro								
Ref.no.	Item		Programs	Temperature	Capacity	Chamber	Dimensions	Price €
<b>8500487</b>	CertoClav Vac Pro 12	Stand-alone internal water tanks	12 programmable	selectable (max. 134°C)	12 l	Stainless steel	20 x 36 cm	<b>3.700,00</b>
<b>8500487W</b>	CertoClav Vac Pro 12	with fixed water connection	12 programmable	selectable (max. 134°C)	12 l	Stainless steel	20 x 36 cm	<b>3.700,00</b>
<b>8500488</b>	CertoClav Vac Pro 18	Stand-alone internal water tanks	12 programmable	selectable (max. 134°C)	18 l	Stainless steel	25 x 35 cm	<b>3.950,00</b>
<b>8500488W</b>	CertoClav Vac Pro 18	with fixed water connection	12 programmable	selectable (max. 134°C)	18 l	Stainless steel	25 x 35 cm	<b>3.950,00</b>
<b>8500489</b>	CertoClav Vac Pro 22	Stand-alone internal water tanks	12 programmable	selectable (max. 134°C)	22 l	Stainless steel	25 x 45 cm	<b>4.200,00</b>
<b>8500489W</b>	CertoClav Vac Pro 22	with fixed water connection	12 programmable	selectable (max. 134°C)	22 l	Stainless steel	25 x 45 cm	<b>4.200,00</b>

Prices EXW excl. VAT

Deliveries are based to the terms and conditions published on our website.

**Accessories for horizontal laboratory autoclaves**

Accessories for CertoClav Vacuum Pro					
Ref.no.	Item			Accessories for	Price €
8500495	Service Set 3		door sealing, bacteria filter, safety valve	for all models	168,00
8500499	CertoClav R/O-System	requires tap water connection	reverse osmosis system for automatic water supply for the autoclave	for models with fixed water connection	590,00
8500510	sediment filter		replace every 6 months	R/O System	25,00
8500511	precarbon filter		replace every 6 months	R/O System	25,00
8500512	membrane filter		replace every 30 months	R/O System	100,00
8500479	tray set with 3 trays	material: Stainless steel	middle tray 18 x 34 cm (W x L) upper and lower tray 14 x 34 cm (W x L)	Vac Pro 12	200,00
8500480	tray set with 3 trays	material: Stainless steel	middle tray 22 x 33 cm (W x L) upper and lower tray 17 x 33 cm (W x L)	Vac Pro 18	260,00
8500481	tray set with 3 trays	material: Stainless steel	middle tray B x L 22 x 43 cm (W x L) upper and lower tray 17 x 43 cm (W x L)	Vac Pro 22	270,00
8500479A	plate for bottles	material: Stainless steel	10 x 34 cm (W x L)	Vac Pro 12	64,00
8500480A	plate for bottles	material: Stainless steel	12 x 33 cm (W x L)	Vac Pro 18	74,00
8500481A	plate for bottles	material: Stainless steel	12 x 43 cm (W x L)	Vac Pro 22	84,00
1300187	log printer		for logging during sterilization on thermal paper	external	309,00
8500496	thermal paper		10 rolls á 40 m	log printer	12,00
8500502	label printer		prints expiration date and protocol information on self-adhesive label	external	298,00
8500507	direct thermal labels	self-adhesive removable	10 rolls	label printer	19,00
8500503	barcode scanner		allows reading the sterilization protocol of an item	external	79,00
8500482	water distiller	fully automatic capacity: 3 l	distills 1.5 liters of water within an hour	external	210,00
8500483	sealing machine	material: Stainless steel max. sealing width: 250 mm	medical packaging sealer for clean and safe sealing of samples and medical instruments sealing temperature adjustable	external	390,00
8500508	CertoServer	local installation of the CertoCloud	communicates with the autoclave via network automatic data backups & browser-based data access	for all models	499,00
8500212	CertoCloud premium license		enables SMS notifications, storage for 100.000 protocols, premium remote support, real-time-data online	for all models	257,50
8500513	IQ/OQ documentation			for all models	800,00
8500514	exhaust filtration	with condensate inactivation	for 200 cycles for safety sterilization of biological hazardous substances	for all models	690,00

Prices EXW excl. VAT

Deliveries are based to the terms and conditions published on our website.