

# CertoClav Price list 2020



Michael Dirix (CEO)  
regional court Linz, FN122912d  
VAT-no.: ATU22821702  
DVR-no.: 0739430, ARA Nr. 950  
EORI no.: ATEOS100004284

Allgemeine Sparkasse OÖ Bank AG  
account number: 0300-003514  
sort code: 20320  
IBAN: AT98 2032 0003 0000 3514  
SWIFT-code: ASPKAT2L

CertoClav Sterilizer GmbH  
Peintner Straße 10  
A-4060 Leonding

+43 732 674 278  
office@certoclav.com  
www.certoclav.com

## Vertical laboratory autoclaves

CertoClav laboratory autoclave EL (Thermostat controlled)							
Ref.no.	Item	Programs	Temperature	Capacity	Chamber	Dimensions	Price €
8501312	<b>CertoClav EL 12 125/140°C</b>	2 fixed	125°C or 140°C	12 l	Aluminum	24 x 26 cm	<b>2.099,00</b>
8501412	<b>CertoClav EL 12 125/140°C – SET</b> with wire basket	2 fixed	125°C or 140°C	12 l	Aluminum	24 x 26 cm	<b>2.212,00</b>
8501312A	<b>CertoClav EL 12 115/121°C</b>	2 fixed	115°C or 121°C	12 l	Aluminum	24 x 26 cm	<b>2.099,00</b>
8501412A	<b>CertoClav EL 12 115/121°C - SET</b> with wire basket	2 fixed	115°C or 121°C	12 l	Aluminum	24 x 26 cm	<b>2.212,00</b>
8501318	<b>CertoClav EL 18 125/140°C</b>	2 fixed	125°C or 140°C	18 l	Aluminum	24 x 39 cm	<b>2.540,00</b>
8501418	<b>CertoClav EL 18 125/140°C - SET</b> with wire basket	2 fixed	125°C or 140°C	18 l	Aluminum	24 x 39 cm	<b>2.616,00</b>
8501318A	<b>CertoClav EL 18 115/121°C</b>	2 fixed	115°C or 121°C	18 l	Aluminum	24 x 39 cm	<b>2.540,00</b>
8501418A	<b>CertoClav EL 18 115/121°C - SET</b> with wire basket	2 fixed	115°C or 121°C	18 l	Aluminum	24 x 39 cm	<b>2.616,00</b>
CertoClav laboratory autoclave MultiControl (Microprocessor controlled with PT100 Sensor)							
Ref.no.	Item	Programs	Temperature	Capacity	Chamber	Dimensions	Price €
8502312A	<b>CertoClav MultiControl 12</b>	6 fixed 2 programmable	selectable (max. 140°C)	12 l	Aluminium	24 x 26 cm	<b>3.689,00</b>
8502412	<b>CertoClav MultiControl 12 - SET</b> with wire basket	6 fixed 2 programmable	selectable (max. 140°C)	12 l	Aluminium	24 x 26 cm	<b>3.754,00</b>
8502318A	<b>CertoClav MultiControl 18</b>	6 fixed 2 programmable	selectable (max. 140°C)	18 l			<b>NOT AVAILABLE ANYMORE</b>
8502418	<b>CertoClav MultiControl 18 - SET</b> with wire basket	6 fixed 2 programmable	selectable (max. 140°C)	18 l			<b>NOT AVAILABLE ANYMORE</b>
Accessories for vertical laboratory autoclaves							
Ref.no.	Item				Accessories for		Price €
8583010	<b>Wire basket with handle</b>	ø 23 cm		Stainless steel 18/10	EL MultiControl		<b>163,90</b>
1250499	<b>Service Set 1</b>	lid gasket, valve gasket, gasket for exhaust cock and manometer, steam valve		gaskets for a lid maintenance	EL MultiControl		<b>45,50</b>
8582020	<b>Condensation vessel</b>	capacity: 5 l		with exhaust hose	EL MultiControl		<b>285,50</b>
8582016	<b>Exhaust hose</b>	length: 1 m			Condensation vessel (8582020)		<b>33,90</b>
8500210	<b>PC software professional version</b>	allows multi-level programs (several time and temperature steps)		unlimited license key	MultiControl		<b>249,50</b>
1300187	<b>Log printer</b>	for logging during sterilization on thermal paper		external device	MultiControl		<b>326,50</b>
8500496	<b>Thermal paper</b>	10 rolls á 40 m			Log printer (1300187)		<b>12,90</b>
8582010	<b>Exhaust filter</b>	with condensation vessel and exhaust hose		0,01 Micron	MultiControl		<b>515,90</b>
8582011	<b>Filter cartridge</b>	0,01 Micron			replacement for exhaust filter (8582010)		<b>272,50</b>

## Vertical laboratory autoclaves

CertoClav laboratory autoclave Classic (Thermostat controlled with timer)							
Ref.no.	Item	Programs	Temperature	Capacity	Chamber	Dimensions	Price €
8510174	<b>CertoClav Classic 125/140°C</b>	2 fixed	125°C or 140°C	18 l	Stainless steel	26 x 35 cm	<b>2.099,00</b>
8510184	<b>CertoClav Classic 115/121°C</b>	2 fixed	115°C or 121°C	18 l	Stainless steel	26 x 35 cm	<b>2.099,00</b>
CertoClav laboratory autoclave MultiControl 2 (freely programmable for any program, Microprocessor controlled with PT100 Sensor)							
Ref.no.	Item	Programs	Temperature	Capacity	Chamber	Dimensions	Price €
8510175	<b>CertoClav MultiControl 2</b> with Controller MC2	6 fixed 2 programmable	selectable (max.140°C)	18 l	Stainless steel	26 x 35 cm	<b>2.999,00</b>
8510298	<b>CertoClav MultiControl 2 Touch</b> with Controller MC2, 7" Touchscreen, Cloud Basic	6 fixed 2 programmable	selectable (max.140°C)	18 l	Stainless steel	26 x 35 cm	<b>3.599,00</b>
CertoClav laboratory autoclave Connect (full automatic laboratory autoclav with cloud connection)							
Ref.no.	Item	Programs	Temperature	Capacity	Chamber	Dimensions	Price €
8510176	<b>CertoClav Connect</b> with 7" Touchscreen, Cloud Basic	6 fixed 2 programmable	selectable (max.140°C)	18 l	Stainless steel	25 x 35 cm	<b>3.993,00</b>
CertoClav laboratory autoclave Cibus (autoclave for sterilization of beverages and food - fast cooling with counter pressure)							
8510185	<b>CertoClav Cibus</b> with 7" Touchscreen, Cloud Basic	6 fixed 2 programmable	selectable (max.140°C)	18 l	Stainless steel	25 x 35 cm	<b>5.999,00</b>
Accessories for vertical laboratory autoclaves							
Ref.no.	Item				Accessories for		Price €
1300149	<b>Wire basket with handle</b>	ø 25 cm		Stainless Steel 18/10	Classic MultiControl 2 MultiControl 2 Touch Connect		<b>163,90</b>
1310186	<b>Service Set 2</b>	lid gasket, control valve gasket, seals for exhaust cock, pressure gauge, steam valve, safety valve		gaskets for a lid maintenance	Classic MultiControl 2 MultiControl 2 Touch Connect		<b>63,50</b>
8582020	<b>Condensation vessel</b>	capacity: 5 l		with exhaust hose	Classic MultiControl 2 MultiControl 2 Touch		<b>285,50</b>
8582016	<b>Exhaust hose</b>	length: 1 m			Condensation vessel (8582020)		<b>33,90</b>
8500210	<b>PC software professional version</b>	allows multi-level programs (several time and temperature steps)		unlimited license key	MultiControl 2		<b>249,50</b>
1300187	<b>Log printer</b>	for logging during sterilization on thermal paper		external device	MultiControl 2 MultiControl 2 Touch		<b>326,50</b>
8500496	<b>Thermal paper</b>	10 rolls á 40 m			Log printer (1300187)		<b>12,90</b>
8582010	<b>Exhaust filter</b>	with condensation vessel and exhaust hose		0,01 Micron	MultiControl 2 MultiControl 2 Touch		<b>515,90</b>
8582011	<b>Filter cartridge</b>	0,01 Micron			replacement for exhaust filter (8582010)		<b>272,50</b>
1300185	<b>Media sensor</b>	for the sterilization of liquids		only available in conjunction with a new device	Connect		<b>374,90</b>
8500508	<b>CertoServer</b>	local server with preinstalled autoclave management software for autoclave management and automatic backups			Connect		<b>823,50</b>
8500212	<b>CertoCloud premium license</b>	enables SMS notifications, storage for protocols, premium remote support, real-time-data online			Connect		<b>269,90</b>

## Accessories for vertical autoclaves

Accessories for vertical autoclaves					
Ref.no.	Item			Accessories for	Price €
8583010	Wire basket with handle	ø 23 cm	Stainless Steel 18/10	Autoclaves aluminum	<b>163,90</b>
1300149	Wire basket with handle	ø 25 cm	Stainless Steel 18/10	Autoclaves S/S	<b>163,90</b>
1252203	Tray	ø 225 mm x 40 mm	Aluminum	Autoclaves aluminum	<b>32,90</b>
1300155	Tray	ø 225 mm x 40 mm	Stainless Steel	Autoclaves S/S	<b>43,90</b>
1250499	Service Set 1	lid gasket, valve gasket, gasket for exhaust cock and manometer, steam valve	gaskets for a lid maintenance	Autoclaves aluminum	<b>45,50</b>
1310186	Service Set 2	lid gasket, control valve gasket, seals for exhaust cock, pressure gauge, steam valve, safety valve	gaskets for a lid maintenance	Autoclaves S/S	<b>63,50</b>
8582020	Condensation vessel	capacity: 5 l	with exhaust hose	all vertical models except Connect	<b>285,50</b>
8582016	Exhaust hose	length: 1 m		Condensation vessel (8582020)	<b>33,90</b>
8500210	PC software professional version	allows multi-level programs (several time and temperature steps)	unlimited license key	MultiControl	<b>249,50</b>
1300187	Log printer	for logging during sterilization on thermal paper	external device	MultiControl MultiControl 2	<b>326,50</b>
8500496	Thermal paper	10 rolls á 40 m		Log printer (1300187)	<b>12,90</b>
8582010	Exhaust filter	with condensation vessel and exhaust hose	0,01 Micron	MultiControl MultiControl 2	<b>515,90</b>
8582011	Filter cartridge	0,01 Micron		replacement for exhaust filter (8582010)	<b>272,50</b>
1300185	Media sensor	for the sterilization of liquids	only available in conjunction with a new device	Connect	<b>374,90</b>
8500508	CertoServer	local server with preinstalled autoclave management software for autoclave management and automatic backups		Connect	<b>823,50</b>
8500212	CertoCloud premium license	enables SMS notifications, storage for protocols, premium remote support, real-time-data online		Connect	<b>269,90</b>
1250000	Check and Service Package			all vertical models	<b>Netto 170,00</b>
1250001	Factory-made validation of the sterilization cycle		with protocol	all vertical models	<b>Netto 549,90</b>
8581000	Modification CertoClav EL	with new controller MC2/PT100		EL	<b>Netto 1.689,50</b>

## Horizontal autoclaves

CertoClav laboratory autoclave Vacuum Pro								
Ref.no.	Item	Programs	Temperature	Capacity	Chamber	Dimensions	Price €	
8500487	<b>CertoClav Vac Pro 12</b>	12 programmable	selectable (max. 134°C)	12 l	Stainless Steel	20 x 36 cm	<b>3.907,00</b>	
8500488	<b>CertoClav Vac Pro 18</b>	12 programmable	selectable (max. 134°C)	18 l	<b>NOT AVAILABLE ANYMORE</b>			
8500489	<b>CertoClav Vac Pro 22</b>	12 programmable	selectable (max. 134°C)	22 l	Stainless Steel	25 x 45 cm	<b>4.434,00</b>	
8500610	<b>CertoClav Vac Pro 45*</b>	<small>*Product set of icanCLAVE autoclave STE-45-T &amp; CertoTouch screen</small>	8 programmable	selectable (max. 134°C)	45 l	Stainless Steel	31,9 x 61,7 cm	<b>9.888,00</b>

## Accessories for horizontal laboratory autoclaves

Accessories for CertoClav Vacuum Pro					
Ref.no.	Item		Material:	Accessories for	Price €
8500479	<b>Tray Set with 3 trays</b>	middle tray 18 x 34 cm (W x L) upper and lower tray 14 x 34 cm (W x L)	Stainless steel	Vac Pro 12	<b>211,50</b>
8500481	<b>Tray Set with 3 trays</b>	middle tray 22 x 43 cm (W x L) upper and lower tray 17 x 43 cm (W x L)	Stainless steel	Vac Pro 22	<b>285,50</b>
8500404	<b>Tray Set with 3 trays</b>	middle tray 22 x 43 cm (W x L) upper and lower tray 17 x 43 cm (W x L)	Stainless steel	Vac Pro 45	<b>325,50</b>
1300765	<b>Drip tray for bottles</b>	298 x 142 x 16 (D x W x H)	Stainless steel	Vac Pro 12	<b>169,90</b>
1300767	<b>Drip tray for bottles</b>	378 x 167 x 16 (D x W x H)	Stainless steel	Vac Pro 22	<b>193,50</b>
1300768	<b>Drip tray for bottles</b>	511 x 222 x 16 (D x W x H)	Stainless steel	Vac Pro 45	<b>210,50</b>
8500401	<b>Wire basket</b>	13,5 x 13,5 x 32 cm (W x H x L)	Stainless steel	Vac Pro 12	<b>122,90</b>
8500403	<b>Wire basket</b>	17 x 17 x 41 cm (W x H x L)	Stainless steel	Vac Pro 22	<b>184,50</b>
8500406	<b>Wire basket</b>	17 x 17 x 41 cm (W x H x L)	Stainless steel	Vac Pro 45	<b>387,50</b>
8500408	<b>SET Wire basket &amp; drip tray</b>		Stainless steel	Vac Pro 12	<b>279,90</b>
8500410	<b>SET Wire basket &amp; drip tray</b>		Stainless steel	Vac Pro 22	<b>369,90</b>
8500411	<b>SET Wire basket &amp; drip tray</b>		Stainless steel	Vac Pro 45	<b>589,90</b>
8500495	<b>Service Set 3</b>	door sealing, bacteria filter, safety valve		Vac Pro 12-22	<b>177,90</b>
8500407	<b>Service Set 5</b>	Door sealing, bacteria filter, safety valve		Vac Pro 45	<b>199,90</b>
1300705	<b>Bacteria filter</b>	replace every 300 cycles		Vac Pro 12-22	<b>56,50</b>
1300187	<b>Log printer</b>	for logging during sterilization on thermal paper	external device	all models	<b>326,50</b>
8500496	<b>Thermal paper</b>	10 rolls á 40 m		Log printer (1300187)	<b>12,90</b>
8500482	<b>Water distiller</b>	distills fully automatic 1.5 liters of water within an hour	external device		<b>221,50</b>
8500502	<b>Label printer</b>	prints expiration date and protocol information on self-adhesive label	external device	all models	<b>314,90</b>

## Accessories for horizontal laboratory autoclaves

8500503	Barcode scanner	allows reading the sterilization protocol of an item	external device	all models	83,50
8500507	Direct thermal labels	10 rolls, self-adhesive, removable		Label printer (8500502)	20,50
8500508	CertoServer	local server with preinstalled autoclave management software for autoclave management and automatic backups		all models	823,50
8500212	CertoCloud premium license	enables SMS notifications, storage for protocols, premium remote support, real-time-data online		all models	269,90
8500513	Template for IQ/OQ documentation	on usb flash drive		Vac Pro 12-22	844,90
8500513	IQ/OQ qualification	Factory-made validation including documentation		Vac Pro 12-22	3.605,00

## Spare parts for CertoClav tabletop autoclaves

Gaskets			
Ref.no.	Item	Spare Part for	Price €
1250404	loose Silicon white 24 cm for lid	Autoclaves aluminum	24,50
1260404	packed Silicon white 24 cm for lid	Autoclaves aluminum	24,90
1300122	white 26 cm for lid	Autoclaves S/S	45,50
1250433	for sterilization valve aluminum lid	Autoclaves aluminum	2,00
1300483	for sterilization valve stainless steel lid	Autoclaves S/S	2,90
1250437	for thermometer for CV-EL	CertoClav EL	4,50
1250451	for manometer and steam release cock	all vertical models	2,90
1250408	for base black 24 cm	Autoclaves aluminum	19,50
1300026	for base white 26 cm	Autoclaves S/S	25,90
Valves			
Ref.no.	Item	Spare Part for	Price €
1250304	CV-EL / Classic 125/140°C	CertoClav EL CertoClav Classic 125/140°C	124,90
1250302	CV-EL / Classic 115/121°C	CertoClav EL CertoClav Classic 115/121°C	124,90
1250320	for MultiControl / Essential 3,0 bar	CertoClav MultiControl CertoClav MultiControl 2 CertoClav Essential	124,90
1250242	Steam release valve aluminum lid	Autoclaves aluminum	31,50
1300123	Steam release valve stainless steel lid	Autoclaves S/S	33,90
1250310	Safety valve complete lid Aluminum	Autoclaves aluminum	31,90
1300103	Safety valve complete lid Stainless steel	Autoclaves S/S	35,00
Handles			
Ref.no.	Item	Spare Part for	Price €
1250235	Lid handle with lock complete	Autoclaves aluminum	45,50
1240239	Locking mechanism 4 parts only for repair	Autoclaves aluminum	32,50
1250236	Lid handle with screw and nut	Autoclaves aluminum	27,90
1250237	Handle with arrow for chamber Aluminum	Autoclaves aluminum	27,50
1250238	Handle without arrow for chamber Alu	Autoclaves aluminum	27,50
1250234	Lid handle with lock complete for ED	Autoclaves S/S	65,90
1300086	Lid handle with screw and nut for ED	Autoclaves S/S	23,90
1300012	Handle with arrow for chamber stainless-steel	Autoclaves S/S	20,50

## Spare parts for CertoClav tabletop autoclaves

Lid			
Ref.no.	Item	Spare Part for	Price €
7940326	Lid Alu without parts CV-EL/MC	CertoClav EL CertoClav MultiControl	429,50
8511310	Lid CV-EL complete, without pressure control valve	CertoClav EL	1.065,90
8510312	Lid MultiControl, ES Alu complete with valve	CertoClav MultiControl CertoClav Essential	1.120,50
Parts			
Ref.no.	Item	Spare Part for	Price €
1250332	Thermometer with gasket	all vertical models	151,50
1250330	Glass scale for thermometer	all vertical models	61,90
1250334A	Manometer EN837-1 Cl.1.6 with Manufact. Certificat acc. EN10204	all vertical models	217,90
1250348	Steam release cock CV-EL complete	CertoClav EL	130,90
1250194	Counter nut for manometer and steam release cock	all vertical models	3,00
1250208	Steam release cock MC with venting nozzle	CertoClav MultiControl CertoClav MultiControl 2	211,90
1240211	Plate with venting nozzle for MC 11 mm ø	CertoClav MultiControl CertoClav MultiControl 2	81,50
1250380	Cleaning needle for venting nozzle	all vertical models	25,90
Aluminum chamber with handles and base			
Ref.no.	Item	Spare Part for	Price €
8520112	CV-EL 12	CertoClav EL 12	1.463,50
8520118	CV-EL 18	CertoClav EL 18	1.985,50
8523312	MultiControl / Essential 12	CertoClav MultiControl 12 CertoClav Essential 12	1.512,00
8523318	MultiControl / Essential 18	CertoClav MultiControl 18 CertoClav Essential 18	2.033,90
Base with socket			
Ref.no.	Item	Spare Part for	Price €
1250324	CV-EL 12/18 with or without fuses	CertoClav EL	130,50
1251324	MultiControl / Essential	CertoClav MultiControl CertoClav Essential	130,50
1300025CL	Classic	CertoClav Classic	147,90
1300025MC	MC2	CertoClav MultiControl 2	147,90



## Spare parts for CertoClav tabletop autoclaves

Cable set complete			
Ref.no.	Item	Spare Part for	Price €
1250357	for CV-EL without fuses	CertoClav EL	59,90
1250356	for MultiControl chamber Aluminum	CertoClav MultiControl	47,50
1310391	for Classic	CertoClav Classic	65,50
1310027	for MC2	CertoClav MultiControl 2	65,50
1250337	Cable with mains plug for CV-EL /Classic 1250337UK = with UK-plug	all vertical models	50,00
Controller			
Ref.no.	Item	Spare Part for	Price €
1310243	Control unit MultiControl 2	CertoClav EL CertoClav MultiControl CertoClav MultiControl 2	1.292,00
1300188K	PT 100 Sensor with cable	CertoClav EL CertoClav MultiControl CertoClav MultiControl 2	233,50
1300197	Mains cable MultiControl 2/Essential	CertoClav EL CertoClav MultiControl CertoClav MultiControl 2	49,90
1300198	Heating cable MultiControl 2/Essential	CertoClav EL CertoClav MultiControl CertoClav MultiControl 2	49,90
1250360	Sensor cable green MC old version	CertoClav MultiControl	192,50
1250381	Mains cable 3 m MultiControl old version	CertoClav MultiControl	49,90
1250382	Heating cable 1,4 m MultiControl old version	CertoClav MultiControl	49,90

## Spare parts for CertoClav tabletop autoclaves

Miscellaneous spare parts			
Ref.no.	Item	Spare Part for	Price €
1250347	Green rocker switch with cover	Autoclaves aluminum	<b>28,90</b>
1250329	Tuning knob setting temp. CV-EL/Classic	CertoClav EL CertoClav Classic	<b>14,00</b>
1250338	Thermostat CV-EL Alu	CertoClav EL	<b>58,90</b>
1230338A	Thermostat Classic Stainless steel	CertoClav Classic	<b>58,90</b>
1250335	Bimetal fuse chamber Aluminum	Autoclaves aluminum	<b>25,50</b>
1300309	Bimetal fuse chamber stainless-steel	Autoclaves S/S	<b>25,50</b>
1250325	Socket complete for base black	Autoclaves aluminum	<b>41,50</b>
1240325A	Socket complete for base white	Autoclaves S/S	<b>41,50</b>
1250363	Fuse 10 Amp.	all vertical models	<b>2,00</b>
1300031	Timer Classic	CertoClav Classic	<b>56,50</b>
1300042	Knob for timer Classic	CertoClav Classic	<b>52,50</b>
1300036	Lamp for Classic	CertoClav Classic	<b>53,50</b>
Repairs			
Ref.no.	Item	Accessories for	Price €
1300000	15 minutes working time	For repairs only original CertoClav spare parts are used. all autoclaves	<b>Netto 33,00</b>



**CERTOCLAV**  
Since 60 years.

**Order**  
**Distributor verification**  
**Custom-made products**  
**Operating instructions and**  
**operating videos**  
**Technical support**



Certoclav Sterilizer GmbH  
Peintner Straße 10  
A-4060 Leonding



Phone: +43 732 674 278  
E-Mail: [support@certoclav.com](mailto:support@certoclav.com)  
Website: [www.certoclav.com](http://www.certoclav.com)